

See information at www.zygem.com
or email support@zygem.com

rapid, enabling solutions
for Biotechnology



DNA Extraction Using *prepGEM™* Insect


Zygem Quick-Start Guide

prepGEM™ Insect

The following method is recommended for extracting DNA from insects or insect parts using *prepGEM™* Insect.

- All manipulations should be performed in a clean-room or a PCR hood.
- Use only certified DNA-free tubes and reagents.
- Wash forceps and scalpels and dissection surfaces in 0.05% Hypochlorite bleach. Rinse thoroughly with DNA-free water.

Procedure Outline



Sample Preparation

- Insects, or parts of insects should be placed in a thin-walled PCR tube for extraction. Approximately 1-2 mm³ of total tissue should be used.
- For older tissue, the sample should be crushed.

To order *prepGEM™* Insect visit www.zygem.com
or email: info@zygem.com or contact your local distributor



Extraction Method

- Add to the material:
35 µl of PCR grade water,
4 µl of 10x Buffer **BLACK**,
1 µl *prepGEM™*

For small insects, these volumes can be reduced further. Solutions can be made as a mastermix and stored on ice.

- Incubate at:
75°C for 15 minutes
95°C for 5 minutes

A thermal cycler can be used for this step

- Transfer supernatant to a new tube

The DNA is in this solution. **Do not discard.**

Centrifugation is undesirable for automation and should not normally be needed. However with some material, two minutes at 13,000 x g may assist in clarifying the extract.

The sample is now ready for PCR.

Depending on the age and quality of tissue, 0.1 - 5 µl of extract is recommended for a 25 µl PCR

Technical tips for sample management

- prepGEM™* Insect is a preparative method for DNA extraction - it is not a purification protocol. The *prepGEM™* method lyses cells and removes nucleoproteins from the DNA. Extracted DNA can be used for many types of genotyping including SNP and STR analysis as well as quantitative, multiplex and end-point PCR.
- DNA extracted using *prepGEM™* is largely single-stranded because of the 95°C heat step.
- For accurate yield assessment, a qPCR is recommended. The DNA produced by *prepGEM™* is approximately 90% single-stranded. If standard fluorescent chelating dyes are to be used for quantification, then this factor should be taken into consideration.
- As with any preparative method for nucleic acid extraction, for best results prepare and manage samples at 4°C, or on ice, before and after extraction.
- Depending on the type of insects present in a sample, and their state of desiccation, grinding and/or re-hydration may be required to achieve optimal yield.
- The presence of EDTA in the sample can reduce enzyme performance in the extraction. This problem can be overcome by adding CaCl₂ to a final concentration of 200 µM.
- When storing sample after extraction, aspirate the supernatant from the tissue and store separately.
- For long term storage of the extracted DNA, add TE buffer to 1x (10 mM Tris, pH 7.5, 1 mM EDTA) and store at -20°C.
- The *prepGEM™* reagents are stable at 4°C but after tubes have been opened and for longer term storage, reagents should be stored at -20°C.
- Zygem Corp Ltd - Research Use Only. All products are subject to a limited use license. See the product documentation and information at www.zygem.com

And then fold here

Fold here