



Antarctica - the source of forensicGEM® (Pressure ridges, Ross Ice Shelf)

PCR inhibition by adhesive tape

Jessica Hayward, University of Auckland, Private Bag 92019, Auckland, New Zealand.

Introduction

Manipulation of trace forensic samples and hair is complicated by the the small size of the sample and static charge generated by the laboratory environment. To simplify collection and handling, samples are often lifted using 3M Post-It® notes or adhesive tape but as a consequence, the tape and sample become inseparable and must be processed together.

The chemical composition of adhesive tape is complex. Each tape contains a different combination of backing film and glue, so it is important to be aware of any of chemical agents that can cause problems for PCR or forensic profiling.

The purpose of the experiment described here, was to assess the level of PCR inhibition caused by several brands of tape when the tape and sample were processed together using the forensicGEM® Tissue one-step DNA extraction kit. To exaggerate the effect of inhibition, extended incubation times were used followed by a PCR inhibition test to obtain a qualitative assessment of inhibition.

Materials and methods

0.25 mm x 0.25 mm sections of adhesive tape (see Figure 1 for tape brands) were cut with a sterile blade and placed into sterile 200 µl PCR tubes. Each sample was washed as follows:

- 200 µl of isopropanol for 5 minutes
- 200 µl of 2.5% Triton X-100 (Merck) for 5 minutes
- 2x washes with DNA/RNA-free water for 5 minutes

The samples were drained and the following extraction reagents added:

1 µl	forensicGEM®
5 µl	10 x Buffer GOLD
44 µl	DNA-free H ₂ O

Samples and reagent blanks were incubated at 75°C for 4 hours and then at 95°C for 15 minutes. Following incubation, samples were centrifuged at 20,000 x g for 5 minutes.

Inhibition PCR

A PCR was performed using AmpliTaq® (Applied Biosystems) in a reaction containing 10 ng of *Escherichia coli* DNA. The primers were designed to amplify the full-length 16S rRNA gene. 10 µl of the tape extract was added to a 50 µl PCR reaction. All reactions were performed in duplicate and the products were visualised on a 0.8% agarose gel.

Results and discussion

The results of the experiment are shown in Figure 1.



Although most forensic DNA extractions do not require prolonged extractions, longer incubations may be necessary for recalcitrant substrates such as bone, teeth and hair shafts. In these situations, it is evident that a careful choice of adhesive tape should be made.

Interestingly, 3M Post-It Notes are a common choice of forensic scientist for tape lifts. The results of this experiment suggest that other brands may be more effective for one-step rapid DNA extraction procedures.

For more information visit: www.zygem.com
or email: info@zygem.com