



## Rapid automated extraction of high quality plant DNA using phytoGEM and the PDQeX<sub>2400</sub>

Dr Abishek Muralidhar

### Introduction

Extracting DNA from plants can be challenging because of their tough cell walls and the inhibitory polysaccharides and tannins released during the extraction process. Currently available methods address these difficulties by using tissue grinding followed by silica beads, columns or organic solvents. These procedures are multistep making them labour intensive and time-consuming.

The *phytoGEM* extraction system from ZyGEM is a simple closed-tube method for extracting PCR-ready plant DNA. It uses a collection tool that crushes the material onto storage cards. Punches from the card are lysed using an aggressive thermophilic proteinase combined with a cocktail of cell wall degrading glycosyl hydrolases. The extraction is performed automatically in the ZyGEM PDQeX<sub>2400</sub> device which both lyses the cells and purifies the DNA through a proprietary polymer matrix.

To demonstrate the extraction efficiency from the *phytoGEM* plant DNA extraction system, six plants were selected: *Arabidopsis thaliana*, *Oryza sativa* (Rice), *Zea mays* (Maize), *Triticum aestivum* (Wheat), *Solanum lycopersicum* (Tomato) and *Citrus limon* (Lemon). Both qPCRs and end-point PCRs were used to evaluate the integrity and quality of the extracted DNA.

### Methods

#### Sample preparation

Healthy leaves from the six plants were crushed onto the storage cards using the *phytoGEM* sampling tool (Fig. 1).

#### DNA extraction

4 mm circular punches were removed from the storage cards using a Harris punch. Two punches were placed into the PDQeX extraction tubes. The following reaction mixture was prepared and dispensed into each tube:

- 2  $\mu$ l *phytoGEM*;
- 2  $\mu$ l *prepGEM*;
- 20  $\mu$ l 5X Buffer;
- 76  $\mu$ l H<sub>2</sub>O

Tubes were sealed by inserting the tapered cap and gently flicked to remove air bubbles and to ensure that the punch discs were completely immersed in the liquid. These were then processed using the PDQeX "Plant" extraction program. The program took a total of 15 minutes.

### Quantification and assessment of the extracted DNA

2  $\mu$ l of extract was used for qPCR amplification in a 10  $\mu$ l reaction containing the following: 5  $\mu$ l Perfecta™ SYBR®Green Fastmix™ (Quanta Biosciences), 0.5  $\mu$ l each of forward and reverse primers and 2  $\mu$ l H<sub>2</sub>O.


Endpoint PCRs were performed to amplify products of different molecular weights from both nuclear and chloroplast genes using 5  $\mu$ l PCR mix (Quanta Biosciences), 0.5  $\mu$ l each of forward and reverse primers and 2  $\mu$ l H<sub>2</sub>O. Amplicons were visualized using agarose gel electrophoresis.

### Results


All *phytoGEM* extractions gave good quality DNA that produced clean qPCR plots (Figure 2). Quantitative Comparison of *phytoGEM* and QIAGEN spin columns showed that *phytoGEM* gave higher DNA yield per mm<sup>2</sup> for most leaf types. (Table 1). The QIAGEN kit, however out performed *phytoGEM* with lemon leaf samples with respect to yield.

Figure 3 shows the average DNA yield for different numbers of punches used in the PDQeX extractions. Yield is proportional to the number of punches. The data demonstrate the flexibility of the method in that it allows the user to decide whether to maximise yield or conserve the sample.

*PhytoGEM*-extracted DNA performs equally well with end-point PCR of larger DNA fragments. Amplicons from nuclear (*adh*, *act2*) and chloroplast genes (*rbcL*, *matK*, *psbA*) of sizes ranging from 200-2300 bp were successfully generated from the six species (Figure 4).



**Figure 1:** The steps for plant DNA extraction using the *phytoGEM* system. Leaf samples are crushed using the *phytoGEM* crusher tool (A) onto storage cards (B). Circular punches are taken (C) and transferred into PDQeX tubes (D). Tubes with punches and reagents are inserted into the PDQeX (E) and run through an extraction program. High quality, clear plant DNA (F) is ready for downstream applications.




**Figure 2.** qPCR plots of alcohol dehydrogenase (*adh*) gene fragments amplified from *phytoGEM* extractions from 6 plants.

| Species                    | Method          | Area (mm <sup>2</sup> ) | Total (ng) | Conc (ng / $\mu$ l) | Recovery (ng/mm <sup>2</sup> ) |
|----------------------------|-----------------|-------------------------|------------|---------------------|--------------------------------|
| Arabidopsis                | QIAGEN          | 1024                    | 1515       | 15.15               | 1.47                           |
| ( <i>A. thaliana</i> )     | <i>phytoGEM</i> | 50.24                   | 94         | 0.94                | 4.34                           |
| Rice                       | QIAGEN          | 958                     | 1791       | 17.91               | 1.87                           |
| ( <i>O. sativa</i> )       | <i>phytoGEM</i> | 50.24                   | 199        | 1.99                | 5.95                           |
| Maize                      | QIAGEN          | 633.94                  | 806        | 8.06                | 1.27                           |
| ( <i>Z. mays</i> )         | <i>phytoGEM</i> | 50.24                   | 120        | 1.2                 | 2.38                           |
| Tomato                     | QIAGEN          | 436.32                  | 1398       | 13.98               | 3.20                           |
| ( <i>S. lycopersicum</i> ) | <i>phytoGEM</i> | 50.24                   | 360        | 3.6                 | 7.17                           |
| Wheat                      | QIAGEN          | 1000                    | 2816       | 28.16               | 2.8                            |
| ( <i>T. aestivum</i> )     | <i>phytoGEM</i> | 50.24                   | 110        | 1.1                 | 2.18                           |
| Lemon                      | QIAGEN          | 679                     | 5019       | 50.19               | 7.39                           |
| ( <i>C. limon</i> )        | <i>phytoGEM</i> | 50.24                   | 123        | 1.23                | 2.45                           |


**Table 1.** Comparison of DNA yields between *phytoGEM* and the Plant Mini Kit from QIAGEN. Extractions were made from 5 different plants. The *phytoGEM* extractions recovered substantially more DNA per mm<sup>2</sup> from *A. thaliana*, *O. sativa*, *S. lycopersicum* and *Z. mays*.

**Discussion**

This study demonstrates the automated *phytoGEM* extraction system and its ability to rapidly produce PCR-ready DNA from both monocots and dicots. The method is fast (~15 minutes) and hands-off making it ideal for high-throughput genetic screening of plants.



**Figure 3.** Concentration of DNA in extracts produced using the PDQeX from different plants using 1 – 4 punches. Data are compiled from 5 extractions from each plant species.



**Figure 4.** Agarose gel images showing amplicons generated from plant extractions. Fragments were generated using primer combinations targeting sequences in 5 genes: *adh* (alcohol dehydrogenase), *act2* (actin 2), *matK* (maturase K), *rbcL* (RuBisCO large subunit) and *psbA* (photosystem II) ranging from ~200–2300 bp.

- (A) *Arabidopsis*: 1-4: *act2*, 5-6: *matK*, 7-8: *rbcL*, 9-10 *psbA*, 11-12: *adh*
- (B) *Rice*: 1-2: *rbcL*, 3-4: *psb*, 5-14: various fragments of *adh*.
- (C) 1-2: lemon *adh*, 3-4: lemon *rbcL*, 5: lemon *psb*, 6-7: maize *adh*, 8-9: maize *psb*, 10-11: maize *rbcL*, 12-13: tomato *psb*, 14-15: tomato *rbcL*.

**References**

- [1] Doyle & Doyle, 1990; *Biotechniques*, 12: 13–15.
- [2] Lickfeldt et. al, 2002; *Hort.Science*, 37: 822–825
- [3] Tamari et. al, 2013; *J Biomol Tech.*, 24: 113–118